

PTGES3 Human

Description: Recombinant Human PTGES3 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 160 amino acids (1-160 a.a.) and having a molecular mass of 18.6 kDa. PTGES3 is purified by conventional chromatography techniques.

Catalog #: ENPS-465

For research use only.

Synonyms: TEBP, CPGES, SID3177, 5730442A20Rik, p23, HSP90 co-chaperone, Prostaglandin E synthase 3, Cytosolic prostaglandin E2 synthase, Telomerase-binding protein p23, Progesterone receptor complex p23, PTGES3.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MQPASAKWYD RRDYVFIEFC VEDSKDVNVN FEKSKLTFSC
LGGSDNFKHL NEIDLFHCID PNDSKHKRTD RSILCCLRKG ESGQSWPRLT KERAKLNWLS
VDFNNWKDWE DDSDEDMNSF DRFSEMMNNM GGDEDVDLPE VDGADDDSDQD
SDDEKMPDLE.

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The PTGES3 protein solution contains 20mM Tris-HCl, pH-8, 1mM DTT & 10% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

PTGES3 takes part as a cochaperone and is involved in signal transduction. PTGES3 is a molecular chaperone that localizes to genomic response elements in a hormone-dependent manner and disrupts receptor-mediated transcriptional activation, by promoting disassembly of transcriptional regulatory complexes. PTGES3 is necessary for appropriate functioning of the glucocorticoid and other steroid receptors. PTGES3 localizes to genomic response elements in a hormone-dependent method and disrupts receptor-mediated transcriptional activation, by promoting disassembly of transcriptional regulatory complexes.

To place an order, please [Click HERE](#).