

MAT2B Human

Description: Recombinant Human MAT2B produced in E.Coli is a single, non-glycosylated polypeptide chain containing 323 amino acids (1-323 a.a) and having a molecular mass of 36.4 kDa. MAT2B is purified by conventional chromatography techniques.

Catalog #: ENPS-468

For research use only.

Synonyms: MAT2-beta, MAT-2B, MAT2-B, DTD-4-keto-6-deoxy-D-glucose 4-reductase, MAT-II, MATIIBeta, MAT II beta, Methionine adenosyltransferase 2 subunit beta, Methionine adenosyltransferase II beta, MGC12237, MSTP045, Nbla02999, SDR23E1, TGR, UNQ2435/PRO4995, MAT2B

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MPMPEDMEQ EEVNIPNRRV LVTGATGLLG RAVHKEFQQN
NWHAVGCGFR RARPKFEQVN LLDSNAVHHI IHDFQPHVIV HCAAERRPDV VENQPDAAASQ
LNVDAAGNLA KEAAVGAFL IYSSDYVFD GTNPPYREED IPAPLNLYGK TKLDGEKAVL
ENNLGAVALR IPILYGEVEK LEESAVTVMF DKVQFSNKSA NMDHWQQRFP THVKDVATVC
RQLAEKRMLD PS

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The MAT2B protein solution contains 20mM Tris-HCl, pH-8, 1mM DTT, 1mM EDTA and 10% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

MAT2B is part of the methionine adenosyltransferase family. MAT2B catalyzes the biosynthesis of S-adenosylmethionine from methionine and ATP. MAT2B is the regulatory beta subunit of MAT. MAT2B expression in hepatoma cell lines increases DNA synthesis and thus take part in cell proliferation.

To place an order, please [Click HERE](#).