

IVD Human

Description: IVD Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 415 amino acids (33-426 a.a.) and having a molecular mass of 45.3 kDa. The IVD is fused to a 20 amino acid his tag at N-terminus and purified by conventional chromatography.

Catalog #: ENPS-497

For research use only.

Synonyms: FLJ12715, isovaleryl-CoA dehydrogenase mitochondrial, FLJ34849, EC 1.3.99.10, IVD, ACAD2.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MHSLLPVDDA INGLSEEQRQ
LRQTMKFLQ EHLAPKAQEI DRSNEFKNLR EFWKQLGNLG VLGITAPVQY GGSGLGYLEH
VLVMEIEISRA SGAVGLSYGA HSNLCINQLV RINGNEAQKEK YLPKLISGEY IGALAMSEPN
AGSDVVSMLK KAEKKGNYHI LNGNKFITN GPDADVLIVY AKTDLAAPVA SRGITAFIVE
KGMPGFSTSK KL

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The IVD protein solution contains 20mM Tris-HCl pH-8, 1mM DTT and 10% glycerol.

Stability:

IVD although stable 4°C for 4 weeks, should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

IVD is a mitochondrial matrix enzyme that is part of the cyl-CoA dehydrogenase family which catalyzes the third step in leucine catabolism. The genetic deficiency of IVD leads to a buildup of isovaleric acid, which is toxic to the central nervous system and results in isovaleric acidemia. IVD is a homotetrameric flavoenzyme which catalyzes the conversion of isovaleryl-CoA to 3-methylcrotonyl-CoA.

To place an order, please [Click HERE](#).