

KLK11 Human

Description: KLK11 Human Recombinant produced in HEK-293 cells is a single, glycosylated, polypeptide chain having a molecular weight of 25.6kDa though on SDS-PAGE migrates at about 41kDa due to the glycosylation. The KLK11 is purified by proprietary chromatographic techniques.

Synonyms: kallikrein-related peptidase 11, TLSP, hippostasin, Kallikrein-11, hK11, Hippostasin, Trypsin-like protease, Serine protease 20, PRSS20, MGC33060, EC 3.4.21.

Source: HEK293-F Cells.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Amino Acid Sequence: ETRIIKGFEC KPHSQPWQAA LFEKTRLLCG ATLIAPRWLL
TAAHCLKPRY IVHLGQHNLQ KEEGCEQTRT ATESFPHPGF NNSLPNKDHR NDIMLVKMAS
PVSITWAVRP LTLSSRCVTA GTSCLSGWG STSSPQLRLP HTLRCANITI IEHQKCNAY
PGNITDTMVC ASVQEGGKDS CQGDSSGGPLV CNQSLQGIIS WQDPCAIR KPGVYTKVCK
YVDWIQETMK NN

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The KLK11 protein was lyophilized from a 0.2

Stability:

Lyophilized KLK11 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution KLK11 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized KLK11 in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Introduction:

KLK11 is a multifunctional protease. KLK11 cleaves "bz-Phe-Arg-4-methylcoumaryl-7-amide", KLK11 is a kallikrein substrate, and weakly cleaves other substrates for kallikrein and trypsin. KLK11 cleaves synthetic peptides after arginine rather than lysine residues.

Biological Activity:

KLK11 specific activity is greater than 2200 pmoles/min/g when measured with 100µM colorimetric peptide substrate (D-Val-Leu-Lys-ThioBenzyl ester) and 1 g of activated enzyme. The reaction is carried out in a volume of 100L containing 50mM Tris, 1M NaCl, 10mM EDTA, 0.1mM DTNB (5,5-Dithio-bis(2-nitrobenzoic acid)), pH 8.5 at 37°C. Cleavage of the peptide substrate can be measured at a wavelength of 405nm and quantified using the extinction coefficient 13,260M⁻¹cm⁻¹.

Catalog #:ENPS-518

For research use only.