

LYPLAL1 Human

Description:LYPLAL1 Human Recombinant produced in E. coli is a single polypeptide chain containing 261 amino acids (1-237) and having a molecular mass of 28.9kDa. LYPLAL1 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:ENPS-651

For research use only.

Synonyms:Lysophospholipase-like protein 1, LYPLAL1.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered clear solution.

Amino Acid Sequence:MGSSHHHHH SSGLVPRGSH MGSMAAASG SVLQRCIVSP
AGRHSASLIF LHGSGDSGQG LRMWIKQLN QDLTFQHIKI IYPTAPPSY TPMKGGISNV
WFDRFKITND CPEHLESIDV MCQVLTDLID EEVKSGIKKN RILIGGFSMG GCMAMHLAYR
NHQDVAGVFA LSSFLNKASA VYQALQKSNG VLPELFQCHG TADELVLHSW AEETNSMLKS
LGVTTKFHSF PN

Purity:Greater than 95% as determined by SDS-PAGE.

Formulation:

The LYPLAL1 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolabs products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Lysophospholipase-like protein 1 (LYPLAL1) is a member of the AB hydrolase 2 family. LYPLAL1 does not exhibit phospholipase nor triacylglycerol lipase activity, able to hydrolyze only short chain substrates as a result of its shallow active site.

To place an order, please [Click HERE](#).