

Cetrorelix

Description: Cetrorelix acetate is a synthetic decapeptide with gonadotropin-releasing hormone (GnRH) antagonistic activity. Cetrorelix acetate is an analog of native GnRH with substitutions of amino acids at positions 1, 2, 3, 6, and 10. The molecular formula is C₇₀H₉₂CIN₁₇O₁₄ (Ac-D-Nal¹-D-Cpa²-D-Pal³-Ser⁴-Tyr⁵-D-Cit⁶-Leu⁷-Arg⁸-Pro⁹-D-Ala¹⁰-NH₂), and the molecular weight is 1431 Dalton, calculated as the anhydrous free base.

Catalog #:HOPS-284

For research use only.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Purity: Greater than 98.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

The Cetrorelix peptide was lyophilized with no additives.

Stability:

Lyophilized Cetrorelix although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Cetrorelix should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Cetrorelix in sterile 18M-cm H₂O not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.

To place an order, please [Click HERE](#).