

Hexarelin

Description: Hexarelin stimulates GH secretion. Hexarelin is more resistant to proteolytic degradation than GHRP-6. Hexarelin has 7 amino acids H-His-D-2-Methyl-Trp-Ala-Trp-D-Phe-Lys-NH₂ and having a molecular weight of 887 Dalton. The Molecular Formula is C₄₇H₅₈N₁₂O₆.

Catalog #:HOPS-295

For research use only.

Physical Appearance: Sterile Filtered lyophilized (freeze-dried) powder.

Purity: Greater than 98.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

The Hexarelin peptide was lyophilized with no additives.

Stability:

Lyophilized Hexarelin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Hexarelin should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Hexarelin 18M-cm H₂O not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.

To place an order, please [Click HERE](#).