

HBsAg adw

Description: HbsAg adw produced Pichia Pastoris, having a molecular weight of approximately 24.0 kDa as shown on SDS-PAGE.

Catalog #: HXPS-879

Source: Pichia Pastoris.

For research use only.

Physical Appearance: Sterile Filtered pale solution.

Purity: Greater than 95.0% as determined by SDS-PAGE Coomassie staining.

Formulation:

Sterile Filtered solution containing 20mM Phosphate Buffer, 154mM sodium chloride.

Stability:

HBsAg Should be stored at 4°C. DO NOT FREEZE.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

To place an order, please [Click HERE](#).