

HBsAg adr

Description: Recombinant HBsAg adr full length monomeric protein contains 227 amino acids of the S-gene a molecular weight of 24 kDa. Small amounts of dimer and trimer forms also exist in the solution.

Catalog #: HXPS-880

Source: Saccharomyces cerevisiae.

For research use only.

Physical Appearance: Sterile Filtered clear solution.

Purity: Greater than 98.0% as determined by analysis by RP-HPLC.

Formulation:

Sterile Filtered solution containing 10mM phosphate, 145mM NaCl and 0.1% NaN₃ as a preservative.

Stability:

HBsAg Should be stored at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

To place an order, please [Click HERE](#).