

HBsAg adw Saccharomyces

Catalog #:HXPS-883

Source:Saccharomyces cerevisae.

Physical Appearance:Sterile Filtered clear solution.

For research use only.

Purity:Greater than 98.0% as determined by SDS-PAGE.

Formulation:

Sterile Filtered solution containing 0.05M phosphate buffer & 0.2M NaCl, pH-7.2.

Stability:

HBsAg protein Should be stored at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

To place an order, please [Click HERE](#).