

HBsAg adw Mutant Q-129-L

Description: Recombinant HBsAg adw antigen was mutated by replacing the Glutamine residue at position 129 with Leucine.

Catalog #: HXPS-889

Source: Pichia pastoris.

For research use only.

Physical Appearance: Sterile Filtered clear solution.

Purity: Greater than 85.0% as determined by SDS-PAGE Coomassie staining.

Formulation:

Sterile Filtered solution containing 20mM Na₂HPO₄, 0.03M NaCl and 0.001% Thimerosal, pH-7.4.

Stability:

HBsAg Should be stored at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

To place an order, please [Click HERE](#).