

GroES Human

Description: Recombinant GroES produced in E.Coli is a single, non-glycosylated polypeptide chain containing 102 amino acids and having a molecular mass of 10 kDa.

Catalog #: HYP5-022

Synonyms: CPN10, GROES, HSP10, HSPE1, Chaperonin-10, 10 kDa heat shock protein mitochondrial, 10 kDa chaperonin, Early-pregnancy factor, EPF.

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MAGQAFRKFL PLFDRVLVER SAAETVTKGG IMLPEKSQ GK
VLQATVVAVG SGSKGKGGEI QPVSVKVGDK VLLPEYGGTK VVLDKDYFL FRDGDILGKY VD.

Purity: Greater than 98.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

The GroES protein contains 20mM Tris buffer pH-8 & 50mM NaCl.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HSP10 is part of the molecular chaperons, that are crucial for their efficient folding of proteins in normal as well as stress conditions. GroES function is to bind to HSP60 in the presence of ATP, thus causing a change in the HSP60 conformation & enclosing the protein substrate within the complex. ATP hydrolysis by chaperonin-60 which destabilizes the HSP10-HSP60 complex, thereby allowing it to dissociate and secrete the substrate protein. GroES having the NCBI accession number of NP_002148 was purified by using conventional chromatography techniques.

To place an order, please [Click HERE](#).