

GroEL Human

Description: Recombinant Human GroEL, HSP60 produced in E.Coli is a single, non-glycosylated polypeptide chain fused to a 20 a.a. His tag at N-terminus containing 593 amino acids (1-573 a.a.) and having a molecular mass of 63kDa. The HSP60 is purified by proprietary chromatographic techniques.

Catalog #: HYPS-023

For research use only.

Synonyms: CPN60, GROEL, HSP60, HSP65, SPG13, CHA60, GROL, crpA, mopA, 60 kDa heat shock protein mitochondrial, Heat shock protein 60, HSP-60, 60 kDa chaperonin, Chaperonin 60, Mitochondrial matrix protein P1, P60 lymphocyte protein, HuCHA60, HSPD1.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MLRLPTVFRQ MRPVSRVLAP
HLTRAYAKDV KFGADARALMLQGVDLLADA VAVTMGPKGR TVIIEQSWGS PKVTKDGVTV
AKSIDLKDKY KNIGAKLVQDVANNTNEEAG DGTATATVLA RSIAKEGFEEK ISKGANPVEI
RRGVMLAVDA VIAELKKQSKPVTTPPEEIAQ VATISANGDK EIGNIISDAM KKVGRKGVIT
VKDGKTLNDE LEIIE

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The GroEL protein contains 20mM Tris-HCl buffer pH-8.0, 5mM DTT and 10% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

GroEL, HSP60 is a chaperonin located in the mitochondria which is responsible for the transportation & refolding of proteins from the cytoplasm directly into the mitochondrial matrix. GroEL is regulated by the HSP10 cochaperonin, which is a single heptameric protein ring having a molecular mass of 10 kDa which form a unique complex with HSP60. HSP10, GroES coordinates the ATPase activity of the HSP60 subunits in order to allow the release of bound polypeptide in a manner that is productive for its correct folding. The GroEL protein having the NCBI accession number NP_002147 was purified by using conventional chromatography techniques.

To place an order, please [Click HERE](#).