

HSPA1B Human

Description: The HSPA1B Protein Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 650 amino acids having a molecular mass of 71.16kDa. HSPA1B is fused to a 10 amino acids His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Synonyms: Heat shock 70kDa protein 1B, HSP70-2, HSP70-1/HSP70-2, FLJ54328, HSPA1A, heat shock 70kDa protein 1A/1B, HSPA1.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Amino Acid Sequence: MKHHHHHHAS AKAAAIGIDL GTTYSCVGVF QHGKVEIIAN
DQGNRTTPSY VAFTDTERLI GDAAKNQVAL NPQNTVFDAK RLIGRKFGDP VVQSDMKHWP
FQVINDGDKP KVQVSYKGET KAFYPEEISS MVLTKMKEIA EAYLGYPVTN AVITVPAYFN
DSQRQATKDA GVIAGLNLVR IINEPTAAAI AYGLDRTGKG ERNVLIFDLG GGTDFVSILT
IDDGIFEVKA TA

Formulation:

The protein was filtered (0.4

Stability:

Lyophilized HSPA1B although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution HSPA1B should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Western blotting, ELISA.

Solubility:

It is recommended to add deionized water to prepare a working stock solution of 0.5mg/ml and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

Introduction:

The HSP70 family is found in many intracellular compartments such as chloroplasts, endoplasmic reticulum, mitochondria, and cytosol. These proteins are induced by a variation of biological stresses, including heat stress, in every organism. HSP70 has a range of functions. For example: HSP70 acts as molecular chaperones facilitating the assembly of multi-protein complexes, HSP70 take part in the translocation of polypeptides across cell membranes and to the nucleus, and HSP70 assists in the proper folding of nascent polypeptide chains. HSP70 is a mitochondrial import equipment and has a crucial part in the cytosolic endoplasmic reticulum. Lately, extracellular localized HSP was discovered to take an important part in the induction of a cellular immune response.

Catalog #:HYP5-028

For research use only.