

H7N7 Netherlands

Description: Recombinant Full-Length H7N7 A/Netherlands/219/03 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells and its Mw is approximately 70,000 Dalton.

Catalog #: IHPS-016

For research use only.

Source: Baculovirus Insect Cells.

Physical Appearance: Sterile Filtered colorless solution.

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The Recombinant H7N7 A/Netherlands/219/03 solution contains 10mM Sodium phosphate, pH 7 and 150mM NaCl.

Stability:

H7N7 A/Netherlands/219/03 Recombinant should be stored at 4°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

To place an order, please [Click HERE](#).