

PDK1 Human

Description: PDK1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 429 amino acids (29-436) and having a molecular mass of 48.6 kDa. PDK1 is fused to 20 a.a. His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PKPS-270

For research use only.

Synonyms: Pyruvate dehydrogenase kinase isoform 1, PDK1, Pyruvate dehydrogenase (lipoamide) kinase isozyme 1 mitochondrial.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MSSDSGSSPA SERGVPGQVD
FYARFSPSPL SMKQFLDFGS VNACEKTSFM FLRQELPVRL ANIMKEISLLPDNLLRTPSV
QLVQSWYIQS LQELDFKDK SAEDAKAIYD FTDTVIRIRN RHNDVIPTMA QGVIEYKESF
GVDPVTSQNV QYFLDRFYMS RISIRMLLNQ HSLLFGGKGK GSPSHRKHIG SINPNCNVLE
VIKIDGYENAR RLC

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

PDK1 solution containing 20mM Tris pH-7, 0.5mM DTT, 0.1M NaCl, 0.1mM EDTA 0.1mM PMSF, 1mM MgCl₂ and 40% glycerol.

Stability:

PDK1 Human Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

PDK1 takes part in the regulation of enzymatic activity of mammalian pyruvate dehydrogenase which is part of a mitochondrial multienzyme complex to catalyze the oxidative decarboxylation of pyruvate and is one of the major enzymes responsible for the regulation of homeostasis of carbohydrate fuels in mammals. PDK1 inhibits glioblastoma growth. PDK1 inhibits the mitochondrial pyruvate dehydrogenase complex by phosphorylation of the e1 alpha subunit, therefore contributing to the regulation of glucose metabolism. PDK1 kinase activity is negatively regulated by binding to 14-3-3.

To place an order, please [Click HERE](#).