

AK4 Human

Description:AK4 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 243 amino acids(1-223a.a.) and having a molecular mass of 27.4kDa. AK4 protein is fused to a 20 amino acid His tag at N-terminus and is purified by standard chromatography.

Catalog #:PKPS-274

For research use only.

Synonyms:Adenylate Kinase 4, AK3, Adenylate Kinase 3-like 1, AK3L1, ATP-AMP transphosphorylase, Adenylate Kinase 3, GTP:AMP phosphotransferase, mitochondrial adenylate kinase-3, AK3L2, MGC166959, EC 2.7.4.3.

Source:Escherichia Coli.

Physical Appearance:Sterile filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MASKLLRAVI LGPPGSGKGT
VCQRIQNFG LQHLSSGHFL RENIKASTEV GEMAKQYIEK SLLVPDHVIT RLMMSELENR
RGQHWLLDGF PRTLGAQAEAL DKICEVDLVI SLNIPFETLK DRLSRRWIHP PSGRVYNLDF
NPPHVHGIDD VTGEPLVQQE DDKPEAVAAR LRQYKDVAKP VIELYKSRGV LHQFSGTETN
KIWPYVYTFLF SN

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

AK4 Human solution (1mg/1ml) containing 20mM Tris HCL pH-8, 2mM DTT, 0.1M NaCl, and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

AK4 belongs to the adenylate kinase family of enzymes. AK4 protein is restricted to the mitochondrial matrix. Adenylate kinases regulate the adenine and guanine nucleotide structures in a cell by catalyzing the reversible transfer of phosphate group among these nucleotides. Hereditary mutations which cause AK4 deficiencies in erythrocytes have been associated to hemolytic anemia.

To place an order, please [Click HERE](#).