

PFKM Human

Description: PFKM Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 800 amino acids (1-780 a.a.) and having a molecular mass of 87.3 kDa. PFKM protein is fused to a 20 amino acid His-Tag at N-terminus and purified by standard chromatography.

Catalog #: PKPS-372

For research use only.

Synonyms: EC 2.7.1.11, GSD7, PFK-1, PFK1, PFKA, PFKX, Phosphofructokinase-M, Phosphofructokinase 1, Phosphohexokinase, Phosphofructo-1-kinase isozyme A, MGC8699, PFKM.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MTHEEHAAK TLGIGKIAIV
LTSGGDAQGM NAAVRVVRV GIFTGARVFF VHEGYQGLVD GGDHIKEATW ESVMMLQLG
GTVIGSARCK DFREREGRLR AAYNLVKRGI TNLCVIGGDG SLTGADTFRS EWSDLLSDLQ
KAGKITDEEA TKSSYLIVGLVGSIDNDFC GTDMTIGTDS ALHRIMEIVD AITTAQSHQ
RTFVLEVMGR HCG

Purity: Greater than 80% as determined by SDS-PAGE.

Formulation:

PFKM Human solution containing 20mM Tris HCL pH-8, 5mM DTT, 0.2M NaCl and 20% glycerol.

Stability:

PFKM human although stable at 4°C for 1 week, should be stored desiccated below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

PFKM is a regulatory glycolytic enzyme that converts fructose 6-phosphate and ATP into fructose 1,6-bisphosphate (through PFK-1), fructose 2,6-bisphosphate (through PFK-2) and ADP. Three phosphofructokinase isozymes exist in humans: muscle, liver and platelet. Mutations in PFKM gene have been related with glycogen storage disease type VII, also identified as Tarui disease.

To place an order, please [Click HERE](#).