

CD83 Human

Description:CD83 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 148 amino acids (20-144a.a) and having a molecular mass of 16kDa. CD83 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-026

For research use only.

Synonyms:BL11, HB15, hCD83, B-cell activation protein, Cell surface protein HB15.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSTPEVKVA CSEDVDLPCT
APWDPQVPYT VSWVKLLEGG EERMETPQED HLRGQHYHQK GQNGSFDAPN ERPYSLKIRN
TTSCNSGTYR CTLQDPDGQR NLSGKVILRV TGCPAQRKEE TFKKYRAE.

Purity:Greater than 90% as determined by SDS-PAGE.

Formulation:

CD83 protein solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl and 30% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

CD83 is a single-pass type I membrane protein and part of the immunoglobulin superfamily of receptors. CD83 involves in the regulation of antigen presentation. A soluble form of CD83 can bind to dendritic cells and inhibit their maturation. There were 3 transcript variants encoding different isoforms for CD83.

To place an order, please [Click HERE](#).