

BUB3 Human

Description:BUB3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 351 amino acids (1-328 a.a.) and having a molecular mass of 39.5kDa.BUB3 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-030

For research use only.

Synonyms:BUB3L, hBUB3, Mitotic checkpoint protein BUB3.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered clear solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSMTGSNEF KLNQPPEDGI
SSVKFSPNTS QFLLVSSWDT SVRLYDVPAN SMRLKYQHTG AVLDCAFYDP THAWSGGLDH
QLKMHDNLNTD QENLVGTHDA PIRCVEYCPE VNVMTGSDW QTVKLWDPRT PCNAGTFSQP
EKVYTLVS VG DRLIVGTAGR RVLVWDLRNM GYVQQRRESS LKYQTRCIRA FPNKQGYVLS
SIEGRVAVEY LD

Purity:Greater than 85.0% as determined by SDS-PAGE.

Formulation:

BUB3 protein solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.4M Urea and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Mitotic checkpoint protein BUB3,(BUB3), is a saved component of the mitotic spindle assembly complex (MCC). The encoded protein has 4 WD repeat domains and has a similar sequence as the yeast BUB3 protein. BUB3 is vital for the kinetochore localization of BUB1 and BUBR1. BUB3 participates in the central spindle checkpoint pathway that operates during early embryogenesis. Furthermore BUB3 has a part in regulating the establishment of correct kinetochore-microtubule attachments.

To place an order, please [Click HERE](#).