

C9ORF103 Human

Description: C9ORF103 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 211 amino acids (1-187) and having a molecular mass of 23.1kDa. C9ORF103 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1032

For research use only.

Synonyms: Chromosome 9 open reading frame 103, bA522I20.2, Gluconate kinase, glucokinase-like protein, probable gluconokinase, GNTK, EC 2.7.1.12.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MGSMAAPGA LLVMGVSGSG
KSTVGALLAS ELGWKFYDAD DYHPEENRRK MGKGIPLNDQ DRIPWLCNLH DILLRDVASG
QRVVLACSAL KKTyrDILTQ GKDGVALKCE ESGKEAKQAE MQLLVVHLSG SFEVISGRLL
KREGHFMPPE LLQSQFETL PPAAPENFIQ ISVDKNVSEI IATIMETLKM K

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

The C9ORF103 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 100mM NaCl, 1mM DTT and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

C9orf103 is a member of the gluconokinase gntK/gntV family. C9orf103 takes part in carbohydrate acid metabolism and D-gluconate degradation.

To place an order, please [Click HERE](#).