

ATOH1 Human

Description: ATOH1 Human Recombinant produced in E. coli is a single polypeptide chain containing 198 amino acids (158-354) and having a molecular mass of 21.3kDa. ATOH1 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-1080

For research use only.

Synonyms: Atonal homolog 1 (Drosophila), Helix-loop-helix protein hATH-1, Class A basic helix-loop-helix protein 14, bHLHa14, HATH1, Math1.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MKQRRLLAANA RERRRMHGLN HAFDQLRNVI PSFNNDKKLS
KYETLQMAQI YINALSELLQ TSPGGEQPPP PPASCKSDHH HLRTAASYEG GAGNATAAGA
QQASGGSQRP TPPGSCRTRF SAPASAGGYS VQLDALHFST FEDSALTAMM AQKNLSPSLP
GSILQPVQEE NSKTSPRSHR SDGEFSPHSH YSDSDEAS

Purity: Greater than 85% as determined by SDS-PAGE.

Formulation:

The ATOH1 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.4M Urea and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

ATOH1 activates E box-dependent transcription in collaboration with TCF3/E47, but the activity is entirely provoked by the negative regulator of neurogenesis HES1. ATOH1 takes part in the differentiation of divisions of neural cells by activating E box-dependent transcription.

To place an order, please [Click HERE](#).