

FUBP1 Human

Description: FUBP1 Human Recombinant produced in E. coli is a single polypeptide chain containing 195 amino acids (279-448) and having a molecular mass of 20.8 kDa. FUBP1 is fused to a 25 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1139

For research use only.

Synonyms: Far upstream element (FUSE) binding protein 1, FUSE-binding protein 1, DNA helicase V, FUBP, FBP, hDH V.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSHMDVPIP RFAVGIVIGR
NGEMIKKIQN DAGVRIQFKP DDGTTPERIA QITGPPDRCQ HAAEIITDLL RSVQAGNPGG
PGPGGRGRGR GQGNWNMGPP GGLQEFNFIV PTGKTGLIIG KGGETIKSIS QQSGARIELQ
RNPPPNADPN MKLFTIRGTP QQIDYARQLI EEKIG

Purity: Greater than 85% as determined by SDS-PAGE.

Formulation:

The FUBP1 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 0.15M NaCl and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

FUBP1 is a ssDNA binding protein which stimulates expression of c-myc in undifferentiated cells and activates the far upstream element (FUSE) of c-myc. Regulation of FUSE by FUBP happens by single-strand binding of FUBP to the non-coding strand. FUBP1 performs as an ATP-dependent DNA helicase.

To place an order, please [Click HERE](#).