

C17ORF49 Human

Description: C17ORF49 Human Recombinant produced in E. coli is a single polypeptide chain containing 195 amino acids (1-172) and having a molecular mass of 20.3 kDa. C17ORF49 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1199

For research use only.

Synonyms: Chromosome 17 open reading frame 49, BPTF-associated protein of 18 kDa, chromatin complexes subunit BAP18, MLL1/MLL complex subunit C17orf49, MGC49942, BAP18.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSMTSASTK VGEIFSAAGA
AFTKLGETLM QLHPVADSSP AGAKWTETEI EMLRAAVKRF GDDLNHISCV IKERTVAQIK
ATVKRKVYED SGIPLPAESP KKGPKKVASG VLSPPPAAPP PSSSSVPEAG GPPIKKQKAD
VTLSALNDS D ANSDVVDIEG LGETPPAKKL NFDQA

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

The C17ORF49 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.1M NaCl, 1mM DTT and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

C17ORF49 is a component of chromatin complexes such as the MLL1/MLL and NURF complexes.

To place an order, please [Click HERE](#).