

Cys-Protein-G His

Description: Cys-Protein G His Tag Recombinant produced in E.Coli is a 202 amino acids protein which contains amino acid 190-384 of the Streptococcus sp with a 6-His tag Cysteine residue at the C-terminus, has a predicted molecular mass of approximately 22.5 kDa. But it migrates with an apparent molecular mass of 32kDa in SDS-PAGE. The Protein G His Tag is purified by proprietary chromatographic techniques.

Catalog #: PRPS-1245

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Purity: Greater than 95.0% as determined by SDS-PAGE and HPLC analyses.

Specificity:

1. Binds with greater affinity to most mammalian immunoglobulins than Protein A, including human IgG3 and rat IgG2a. 2. Does not bind to human IgM, IgD and IgA.

Formulation:

Lyophilized with no additive.

Stability:

Lyophilized Protein G although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Protein G should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Protein G binds to the constant region of many species of immunoglobulin G. It can be used to detect, quantify and purify IgG antibodies and antibody/antigen complexes. Recombinant Protein G contains only IgG binding domains. The albumin-binding domain as well as cell wall and cell membrane binding domains have been removed to ensure the maximum specific IgG binding capacity.

Solubility:

It is recommended to reconstitute the lyophilized Protein G in sterile 18M-cm H₂O not less than 100

To place an order, please [Click HERE](#).