

ASGR2 Human

Description:ASGR2 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 255 amino acids (80-311 a.a) and having a molecular mass of 28.9kDa.ASGR2 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-1268

For research use only.

Synonyms:Asialoglycoprotein receptor 2, ASGP-R 2, ASGPR 2, C-type lectin domain family 4 member H2, Hepatic lectin H2, HL-2, ASGR2, CLEC4H2, HBXBP.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSQSEGHRG AQLQAE LRSL
KEAFSNFSSS TLTEVQAIST HGGSVGDKIT SLGAKLEKQQ QDLKADHDAL LFHLKHFPVD
LRFVACQMEL LHSNGSQRTC CPVNWVEHQG SCYWFSHSGK AWAEAEKYCQ LENAHLVVIN
SWEEQKFIVQ HTNPFNTWIG LTDSDGSWKW VDGTDYRHNY KNWAVTQPDN WHGHELGGSE
DCVEVQPDGR WN

Purity:Greater than 90.0% as determined by SDS-PAGE.

Formulation:

ASGR2 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl, 30% glycerol and 1mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Asialoglycoprotein Receptor 2 (ASGR2) is a hetero-oligomeric protein comprised of major and minor subunits, encoded by different genes. ASGR2 is a subunit of the asialoglycoprotein receptor; it is the less abundant minor subunit. The ASGR2 receptor is a transmembrane protein, which has a crucial role in serum glycoprotein homeostasis by mediating the endocytosis and lysosomal degradation of glycoproteins with exposed terminal galactose or N-acetylgalactosamine residues.

To place an order, please [Click HERE](#).