

FIMC E.Coli

Description: FIMC E.Coli Recombinant fused with a 21 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 226 amino acids (37-241 a.a.) and having a molecular mass of 25kDa. The FIMC is purified by proprietary chromatographic techniques.

Catalog #: PRPS-127

For research use only.

Synonyms: Chaperone protein fimC, fimC, b4316, JW4279.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGVALGATRV IYPAGQKQEQ
LAVTNNDENS TYLIQSWVEN ADGVKDGRFI VTPPLFAMKG KKENTLRILD ATNNQLPQDR
ESLFWMNVKA IPSMDKSKLT ENTLQLAIIS RIKLYYRPAK LALPPDQAAE KLRFRRSANS
LTLINPTPY LTVTELNAGT RVLENALVPP MGESTVKLPS DAGSNITYRT INDYGALTPK
MTGVME.

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

The FIMC solution (1 mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 10% glycerol and 1mM DTT.

Stability:

FIMC should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

The 23kDa two-domain periplasmic chaperone FIMC from E.coli is vital for the assembly of type-1 pili, which are filamentous, very oligomeric protein complexes anchored to the outer bacterial membrane that mediate adhesion of pathogenic E. coli strains to host cell surfaces and persist in macrophages.

To place an order, please [Click HERE](#).