

KIN Human

Description: KIN Human Recombinant produced in E. coli is a single polypeptide chain containing 416 amino acids (1-393) and having a molecular mass of 47.8 kDa. KIN is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1282

For research use only.

Synonyms: BTCD, KIN17, KIN, Binding to curved DNA, antigenic determinant of recA protein homolog, DNA/RNA-binding protein KIN17.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MGSMGKSDFL TPKAIANRIK
SKGLQKLRWY CQMCQKQCRDENGFKCHCMS ESHQRQLLLA SENPQQFMDY FSEEFNRNDFL
ELLRRRFGTK RVHNNIVYNE YISHREHIHM NATQWETLTD FTKWLGREGL CKVDETPKGW
YIQYIDRDPE TIRRQLELEKKKKQDLDEE KTAKFIEEQV RRGLEGKEQE VPTFTELSRE
NDEEKVTFNL SKGA

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The KIN solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

DNA/RNA-binding protein KIN17 (KIN17) is a nuclear protein which constructs intranuclear foci for the duration of proliferation and is reallocated in the nucleoplasm during the cell cycle. KIN17 is ubiquitously expressed having the highest levels of expression in the muscle, heart and testis. SV40- altered fibroblasts overexpress KIN17, which interacts with Large T antigen and decreases T-antigen-dependent DNA replication.

To place an order, please [Click HERE](#).