

TNNI2 Human, His

Description: TNNI2 Human Recombinant produced in E. coli is a single polypeptide chain containing 202 amino acids (1-182) and having a molecular mass of 23.5 kDa. TNNI2 is fused to 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1290

For research use only.

Synonyms: AMCD2B, DA2B, FSSV, fsTnI, Troponin I, fast skeletal muscle, Troponin I, fast-twitch isoform.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGDEEKRNR ITARRQHLKS
VMLQIAATEL EKEESRREAE KQNYLAEHCP PLHIPGSMSE VQELCKQLHA KIDAAEEEEKY
DMEVRVQKTS KELEDMNQKL FDLRGKFKRP PLRRVRMSAD AMLKALLGSK HKVCMDLRAN
LKQVKKEDTE KERDLRDVGD WRKNIEEKSG MEGRKKMFES ES.

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

The TNNI2 solution (0.25mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 2M Urea, 20% glycerol and 0.2M NaCl.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Troponin I, fast skeletal muscle (TNNI2), belongs to the troponin I gene family, and a part of the troponin complex including troponin T, troponin C and troponin I subunits. The troponin complex, along with tropomyosin, is responsible for the calcium-dependent regulation of striated muscle contraction. TNNI2 is also suppresses tumor growth in human ovarian carcinoma. Mutations in TNNI2 cause myopathy and distal arthrogryposis type 2B.

To place an order, please [Click HERE](#).