

Protein A

Description: Recombinant Staphylococcal Protein A produced in E.Coli is a non-glycosylated, Polypeptide chain having a molecular mass of 45 kDa. Recombinant Staphylococcal Protein A is purified by proprietary chromatographic techniques.

Catalog #: PRPS-363

Synonyms: Immunoglobulin G-binding protein A, IgG-binding protein A, Staphylococcal protein A, SPA.

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear colorless solution.

Purity: Greater than 98.0% as determined by RP-HPLC.

Formulation:

The protein solution contains no additives.

Stability:

SPA should be stored at -20°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Protein A is a cell wall protein deriving from Staphylococcus aureus which exhibits unique binding properties for IgG from a variety of mammalian species and for some IgM and IgA as well. It binds with the Fc region of immunoglobulins through interaction with the heavy chain. It couples to a wide variety of reporter molecules including fluorescent dyes, enzyme markers, biotin, colloidal gold and radioactive iodine without affecting the antibody binding site. Recombinant Protein A was developed to increase the specificity of the molecule for IgG and is widely used both in research and bioprocessing. The recombinant protein A is produced by expressing a modified protein A gene in E.coli. A specific purification process with strict quality control was taken to get the recombinant protein A with the purity of more than 98% , no human IgG affinity step is used during validated fermentation and purification and devoid of bacterial contaminant found normally in native Protein A. (Free of Staphylococcus endotoxins and hemolysin).

Biological Activity:

Greater than 95.0% binding activity to human IgG.

To place an order, please [Click HERE](#).