

## SNCB Human

**Description:** Recombinant Human b-Synuclein produced in E.Coli is a single, non-glycosylated polypeptide chain (1-134 a.a) having a molecular mass of 14.2kDa (molecular size on SDS-PAGE will appear higher). The Recombinant Human b-Synuclein is purified by proprietary chromatographic techniques.

Catalog #: PRPS-401

For research use only.

**Synonyms:** Beta-synuclein, SNCB, b-Synuclein.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile filtered colorless solution.

**Amino Acid Sequence:** MDVFMKGLSM AKEGVVAAAE KTKQGVTEAA EKTKEGVLYV  
GSKTREGVVQ GVASVAEKT EQASHLGGAV FSGAGNIAAA TGLVKREEFP TDLKPEEVAQ  
EAAEEPLIEP LMEPEGESYE DPPQEEYQEY EPEA.

**Purity:** Greater than 95.0% as determined by SDS-PAGE.

**Formulation:**

The protein (1mg/ml) contains 20mM Tris-HCl buffer (pH 7.5), 0.1M NaCl and 1mM MgCl<sub>2</sub>.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

b-synuclein (amino acids 1-134), an acidic neuronal protein of 134 amino acids, is extremely heat-resistant. b-synuclein also has a chaperone-like activity. Recently, b-synuclein has been suggested to inhibit a-synuclein aggregation. b and a-synuclein are abundantly expressed in the brain and putatively inhibit phospholipase D2 selectively. The encoded protein, which may play a role in neuronal plasticity, is abundant in neurofibrillary lesions of patients with Alzheimer disease.

**To place an order, please [Click HERE](#).**