

Avidin Biotin

Description: Avidin is a basic charged glycoprotein present in the egg white and in some extent in tissues of various animals. Avidin is a homotetrameric protein that contains 4 identical subunits which bind to Biotin with high extent of affinity and specificity consisting of 128 amino acids and has a molecular mass of about 67K Daltons. Avidin-Biotin is a modified form of affinity purified avidin that combines the high affinity specific activity of the native egg white avidin to biotin molecule together with a very low background of the streptavidin produced by the bacteria *Streptomyces avidinii*.

Catalog #: PRPS-453

For research use only.

Synonyms: Avidin, AVD, AVID.

Physical Appearance: Lyophilized freeze dry Powder.

Formulation:

Salt-free lyophilized powder with no preservative added.

Stability:

Store protein solution at 2-8°C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

We highly recommended to prepare a 1 mg/ml solution; Reconstitute using cold deionized pyrogen free water and mix gently until powder is completely dissolved.

Biological Activity:

Binding activity of Avidin-Biotin is 10-15 units/mg protein. One unit of Avidin activity is defined as the amount of protein which will bind 1

To place an order, please [Click HERE](#).