

Bcl 6 Human

Description: Bcl-6 Human Recombinant (aa 1-150) expressed in E.coli, shows a 43kDa band on SDS-PAGE (Including GST tag). The Bcl-6 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-467

Synonyms: B-cell lymphoma 6 protein, BCL-6, Zinc finger protein 51, LAZ-3 protein, BCL-5, Zinc finger and BTB domain-containing protein 27, BCL6, BCL5, LAZ3, ZBTB27, ZNF51, BCL6A.

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Formulation:

Bcl-6 protein in 50mM Tris-Acetate, pH7.5, 1mM EDTA, 20% Glycerol.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

ELISA Inhibition Assays Western Blotting.

Introduction:

Bcl-6 proto-oncogene product (a Kruppel-type zinc-finger protein) is mainly expressed in normal germinal center B cells and related lymphomas. Bcl-6 is involved in chromosome rearrangements at 3q27 in non-Hodgkins lymphomas and Bcl-6 rearrangements have been detected in 33%-45% of diffuse large B cell lymphomas. Bcl-6 has been detected immunohistochemically in follicular lymphomas, diffuse large B cell lymphomas, Burkitts lymphomas and in nodular, lymphocyte predominant Hodgkins disease.

To place an order, please [Click HERE](#).