

TPM1 Human

Description: TPM1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 304 amino acids (1-284 a.a.) and having a total molecular mass of 35kDa (Molecular weight on SDS-PAGE will appear higher). TPM1 is fused to a 20 amino acid His Tag at N-terminus and is purified by proprietary chromatographic techniques.

Synonyms: Tropomyosin alpha-1 chain, Tropomyosin-1, Alpha-tropomyosin, TPM1, C15orf13, TMSA, CMD1Y, HTM-alpha.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MDAIKKKMQM LKLDKENALD
RAEQAEADKK AAEDRSKQLE DELVSLQKKL KGTEDELDKY SEALKDAQEKLELAEEKKATD
AEADVASLNR RIQLVEEELD RAQERLATAL QKLEEAKEAA DESERGMKVI ESRAQKDEEK
MEIQEIQLKE AKHIAEDADRKYEEVARKLV IIESDLERAE ERAELSEGQV RQLEEQLRIM
DQTLKALMAA EDKY

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

TPM1 0.5mg/ml protein solution contains 20mM Tris-HCl buffer pH-8, 1mM DTT, 0.1M NaCl & 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

TPM1 is a member of the tropomyosin family which consists of a number of extremely conserved, extensively distributed 35-45 kDa actin-binding proteins that are involved in the contractile system of striated and smooth muscles and the cytoskeleton of non-muscle cells. Tropomyosin-1 is composed of 2 alpha-helical chains arranged as a coiled-coil. TPM1 is polymerized end to end alongside the two grooves of actin filaments and provides stability to the filaments. TPM1 binds to actin filaments in muscle and non-muscle cells. TPM1 also functions in association with the troponin complex to regulate the calcium-dependent interaction of actin and myosin during muscle contraction. In non-muscle cells TPM1 is implicated in stabilizing cytoskeleton actin filaments. Smooth muscle contraction is controlled by interaction with caldesmon. Alternatively spliced transcript variants encoding a range of isoforms have been described in smooth muscle and non-muscle cells. TPM1 Isoform 1 is expressed in adult and fetal skeletal muscle and cardiac tissues, with higher expression levels in the cardiac tissues, whereas Isoform 10 is expressed in adult and fetal cardiac tissues, but not in skeletal muscle. Mutations in the TPM1 gene are linked to type 3 familial hypertrophic cardiomyopathy.

Catalog #:PRPS-476

For research use only.