

## DYNLRB1 Human

**Description:**DYNLRB1 Recombinant Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 104 amino acids (1-96 a.a.) and having a molecular mass of 11.9 kDa. The DYNLRB1 is fused to 8 amino acid His-Tag at C-terminus and purified by proprietary chromatographic techniques.

**Catalog #:**PRPS-496

For research use only.

**Synonyms:**BITH, BLP, DNCL2A, DNLC2A, ROBLD1, HSPC162.

**Source:**Escherichia Coli.

**Physical Appearance:**Sterile filtered colorless solution.

**Amino Acid Sequence:**MAEEVEETLKR LQSQKGVQGI IVVNTGPII KSTMDNPPTT  
QYASLMHSFI LKARSTVRDIDPQNDLTFLR IRSKKNEIMV APDKDYFLIVIQNPTELEHH HHHH.

**Purity:**Greater than 95.0% as determined by SDS-PAGE.

**Formulation:**

1mg/ml solution containing 20mM Tris pH-8, & 10% glycerol.

**Stability:**

DYNLRB1 Human although stable at 4C for 1 week, should be stored below -18C. Please prevent freeze thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

Dynein light chain roadblock-type (DYNLRB1) belongs to the roadblock dynein light chain family and encodes a cytoplasmic protein which is capable of binding intermediate chain proteins. Upregulation of the DYNLRB1 gene is linked with hepatocellular carcinomas, suggesting that DYNLRB1 may be involved in tumor progression.

**To place an order, please [Click HERE](#).**