

CRMP-1 Human

Description: CRMP-1 Human Recombinant (aa 448-568) expressed in E.coli, shows a 39 kDa band on SDS-PAGE (Including GST tag). The CRMP-1 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-500

For research use only.

Synonyms: Dihydropyrimidinase-related protein 1, DRP-1, Collapsin response mediator protein 1, CRMP-1, CRMP1, DPYSL1, DRP1.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Formulation:

CRMP-1 in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

ELISA Inhibition Assays Western Blotting.

Introduction:

ion and demonstrates that CRMP1 is a bilobed lung-shaped protein forming a tetrameric assembly.

To place an order, please [Click HERE](#).