

FSCN1 Human

Description: FSCN1 Recombinant Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 513 amino acids (1-493 a.a.) and having a molecular mass of 56.7 kDa. The FSCN1 is fused to a 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PRPS-534

For research use only.

Synonyms: Strongylocentrotus Purpuratus, actin-bundling protein, FAN1, HSN, SNL.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MTANGTAEAV QIQFGLINCG
NKYLTAEEAFG FKVNASASSL KKKQIWTLEQ PPDEAGSAAV CLRSHLGRL AADKDGNTVC
EREVPGPDCR FLIVAHDDGR WSLQSEAHRR YFGGTEDRLS CFAQTVSPAE KWSVHIAMHP
QVNIYSVTRK RYAHLSARPA DEIAVDRDVP WGVDSLITLA FQDQRYSVQT ADHRFLRHDG
RLVARPEPAT GY

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

1mg/ml solution containing 20mM Tris pH-8, 2mM DTT, 0.1M NaCl & 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

FSCN1 is an actin-bundling protein that is responsible for rigidity to filopodial bundles to efficiently push the membrane forward during cytoskeleton remodeling and cell migration. FSCN1 takes part in the gathering of actin filament bundles present in microspikes, membrane ruffles, and stress fibers. FSCN1 is absent from most normal epithelia though is upregulated in multiple forms of human carcinoma, where its expression associates clinically with a poor prognosis. FSCN1 expression is abundantly found in neurons, glial cells, and endothelial cells.

To place an order, please [Click HERE](#).