

CBFB HUman

Description: CBFB Recombinant Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 182 amino acids (1-202 a.a.) and having a molecular mass of 23.6 kDa. The CBFB is fused to 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PRPS-541

For research use only.

Synonyms: PEBP2B, polyomavirus enhancer binding protein b, PEA2, CBF-beta.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MPRVVPDQRS KFENEFFRK
LSRECEIKYT GFRDRPHEER QARFQNACRD GRSEIAFVAT GTNLSLQFFPASWQGEQRQT
PSREYVDLER EAGKVYLKAP MILNGVCVIW KGWIDLQRLD GMGCLEFDEE RAQQEDALAQ
QAFEEARRRT REFEDRDRSHREEMEVRSQ LLAVTGKKT RP.

Purity: Greater than 85.0% as determined by SDS-PAGE.

Formulation:

0.5mg/ml solution containing 20mM MES pH-6, 0.1mM PMSF & 10% glycerol.

Stability:

CBFB Human although stable at 4C for 1 week, should be stored below -18C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

CBFB beta subunit is a heterodimeric core-binding transcription factor that is part of the PEBP2/CBF transcription factor family which controls a host of genes particularly to hematopoiesis and osteogenesis. CBFB is a non-DNA binding regulatory subunit which increases DNA binding by alpha subunit.

To place an order, please [Click HERE](#).