

Calumenin Human

Description: Recombinant Human Calumenin produced in E.Coli is a single, non-glycosylated polypeptide chain containing 317 amino acids (20-315 a.a) and having a molecular mass of 37.2 kDa. Calumenin is fused to a 20 amino acid His Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PRPS-703

For research use only.

Synonyms: CALU, Crocalbin, IEF SSP 9302, FLJ90608, Calumenin.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MKPTEKKDRV HHEPQLSDKV
HNDASQFDYD HDAFLGAEAA KTFDQLTPEE SKERLGKIVS KIDGDKDGFV TVDELKDWIK
FAQKRWIYED VERQWKGHDL NEDGLVSWEE YKNATYGYVL DDPDPDDGFN YKQMMVRDER
RFKMADKDG D LIATKEEFTA FLHPPEYDYM KDIVVQETME DIDKNADGFI DLEEYIGDMY
SHDGNTDEPE WV

Purity: Greater than 90% as determined by Analysis by SDS-PAGE.

Formulation:

The Calumenin protein contains 20mM Tris-HCl buffer pH-8 and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Calumenin is a calcium-binding protein located in the endoplasmic reticulum (ER) and sarcoplasmic reticulum (SR) of mammalian tissues which plays a role in ER functions as protein folding and sorting. Calumenin belongs to a family of multiple EF-hand proteins (CERC) that include reticulocalbin, ERC-55, and Cab45 and the product of this gene. Calumenin binds 7 calcium ions having low affinity and takes part in such ER functions as protein folding and sorting.

To place an order, please [Click HERE](#).