

HINT1 Human

Description: HINT1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 126 amino acids (1-126 a.a.) and having a molecular mass of 13.8 kDa. The HINT1 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-709

For research use only.

Synonyms: HINT, PKCI-1, PRKCNH1, FLJ30414, FLJ32340, HINT1, Histidine triad nucleotide-binding protein 1, Adenosine 5'-monophosphoramidase, Protein kinase C inhibitor 1, Protein kinase C-interacting protein 1, PKCI1.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MADEIAKAQV ARPGGDTIFG KIIRKEIPAK IIFEDDRCLA FHDISPQAPT
HFLVIPKKHI SQISVAEDDD ESSLGHLMIK GKKCAADLGL NKGYRMVVNE GSDGGQSVYH
VHLHVLGGRQ MHWPPG.

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The HINT1 solution contains 20mM Tris pH-8 & 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

HINT1, also known as Histidine triad nucleotide-binding protein 1 is part of the superfamily named for a near C-terminal HXHXHXX motif (H: Histidine, X: a hydrophobic amino acid) positioned at the -phosphate of nucleotide substrates. HINT1 hydrolyzes adenosine 5'-monophosphoramidate substrates such as AMP-morpholidate, AMP-N-alanine methyl ester, AMP-alpha-acetyl lysine methyl ester and AMP-NH₂. Though it was initially considered to be a protein kinase C inhibitor and act as a haplo-insufficient tumor suppressor including spontaneous tumor formation in Hint^{+/+} and Hint^{-/-}, its actual physiologic function is not known.

To place an order, please [Click HERE](#).