

CASQ2 Human

Description: CASQ2 Human Recombinant fused to 37 amino acid His Tag at N-terminal produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 417 amino acids (20-399 a.a.) and having a molecular mass of 48.4 kDa. The CASQ2 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-806

For research use only.

Synonyms: PDIB2, CASQ2, Calsequestrin-2, Calsequestrin cardiac muscle isoform, FLJ26321, FLJ93514.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear colorless solution.

Amino Acid Sequence: MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMEEG
LNFTPTDYGKD RVVSLSEKNF KQVLKDYDLL CLYYHEPVSS DKVTQKQFQL KEIVLELVAQ
VLEHKAIGFV MVDAAKEAKL AKKLGDFDEEG SLYILKGDRIT IEFDFGEFAAD VLVEFLDLI
EDPVEIISSEK LEVQAFERIE DYIKLIGFFK SEDSEYYKAF EAAAEHFQPY IKFFATFDKG
VAKKLSLKMN EV

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

The CASQ2 solution contains 20mM Tris-HCl pH-8, 0.1M NaCl, 1mM DTT and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please avoid freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

CASQ2 belongs to the calsequestrin family and is localized to the sarcoplasmic reticulum in cardiac and slow skeletal muscle cells. CASQ2 is a calcium binding protein that stores calcium for muscle function. The discharge of calcium bound to CASQ2 through a calcium release channel activates muscle contraction. CASQ2 binds 40 to 50 moles of calcium. CASQ2 mutations result in stress-induced polymorphic ventricular tachycardia, also called catecholaminergic polymorphic ventricular tachycardia 2 which is known for its bidirectional ventricular tachycardia that causes cardiac arrest.

To place an order, please [Click HERE](#).