

RAB13 Human

Description: RAB13 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 220 amino acids (1-200) and having a molecular mass of 24.7 kDa. The RAB13 is fused to a 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: PRPS-889

For research use only.

Synonyms: RAB13 member RAS oncogene family, Cell growth-inhibiting gene 4 protein, Growth-inhibiting gene 4 protein, RAS-associated protein RAB13.

Source: Escherichia Coli.

Physical Appearance: RAB13 is supplied as a sterile filtered clear solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MAKAYDHLFK LLLIGDSGVG
KTCLIRFAE DNFNNTYIST IGIDFKIRTV DIEGKKIKLQ VWDTAGQERF KTITTAYYRG
AMGIILVYDI TDEKSFENIQ NWMKSIKENA SAGVERLLLG NKCDMEAKRK VQKEQADKLA
REHGIRFFET SAKSSMNVDE AFSSLARDIL LKSGGRRSGN GNKPPSTD LK TCDKKNTNKC

Purity: Greater than 85.0% as determined by SDS-PAGE.

Formulation:

RAB13 protein (0.25mg/ml) is supplied in 20mM Tris-HCL, pH-8, 1mM DTT, 200mM NaCl, 1mM EDTA and 20% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

RAB13 is a member of the small GTPase superfamily which are nontransforming monomeric GTP-binding proteins that hold 4 extremely conserved regions that take part in GTP binding and hydrolysis. In addition, RAB proteins have a vital role in the endocytic pathways. The GTPase Rab13 regulates the assembly of functional epithelial tight junctions (TJs) and is expressed in every phase of preimplantation development.

To place an order, please [Click HERE](#).