

SEPT6 Human

Description: SEPT6 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 457 amino acids (1-434) and having a molecular mass of 52.1 kDa. SEPT6 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-956

For research use only.

Synonyms: Septin 6, Septin 2, SEP2, KIAA0128, SEPT2.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSMAATDIA RQVGEGCRTV
PLAGHVGFDS LPDQLVNKSV SQGFCFNILC VGETGLGKST LMDTLFNTKF EGEPATHTQP
GVQLQSNTYD LQESNVRLKL TIVSTVGFGD QINKEDSYKP IVEFIDAQFE AYLQEELKIR
RVLHTYHDSR IHVCLYFIAP TGHSLKSLDL VTMKKLDSKV NIPIIAKAD AISKSELTKE
KIKITSELVS NG

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

The SEPT6 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 0.1M NaCl and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

SEPT6 is a member of the septin family of GTPases and is functionally involved in cytokinesis and conservation of cellular morphology. SEPT6 interact with SEPT2. The genes encoding the mixed-lineage leukemia and SEPT2 breakpoint proteins are associated with one type of pediatric acute myeloid leukemia which is caused by reciprocal translocation between chromosomes 11 and X.

To place an order, please [Click HERE](#).