

## ARL4D Human

**Description:**ARL4D Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 221 amino acids (1-201) and having a molecular mass of 24.3 kDa.ARL4D is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

**Catalog #:**PRPS-960

For research use only.

**Synonyms:**ADP-ribosylation factor-like 4D, ADP-ribosylation factor-like protein 4L, ARF4L, ARL6.

**Source:**E.coli.

**Physical Appearance:**Sterile Filtered colorless solution.

**Amino Acid Sequence:**MGSSHHHHHH SSGLVPRGSH MGNHITEMAP TASSFLPHFQ  
ALHVVVIGLD SAGKTSLLYR LKFKEFVQSV PTKGFNTEKI RVPLGGSRGI TFQVWDVGGQ  
EKLRPLWRSY TRRTDGLVFV VDAAEAERLE EAKVELHRIS RASDNQGVVPV LVLANKQDQP  
GALSAAVEK RLAVRELAAL TLTHVQGCSA VDGLGLQQGL ERLYEMILKR KKAARGGKKR R

**Purity:**Greater than 90% as determined by SDS-PAGE.

### Formulation:

The ARL4D solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 0.2M NaCl and 20% glycerol.

### Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Introduction:

ARL4D is a developmentally regulated member of the ADP-ribosylation factor/ARF-like protein (ARF/ARL) family of Ras-related GTPases. Small GTP-binding protein which cycles between an inactive GDP-bound and an active GTP-bound form, and the rate of cycling is regulated by guanine nucleotide exchange factors (GEF) and GTPase-activating proteins (GAP). ARL4D takes part in membrane-associated intracellular trafficking. Mutations in this gene is linked to Bardet-Biedl syndrome (BBS).

**To place an order, please [Click HERE](#).**