

ARL15 Human

Description:ARL15 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 224 amino acids (1-204) and having a molecular mass of 25.0 kDa. ARL15 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-963

For research use only.

Synonyms:ADP-ribosylation factor-like 15, ADP-ribosylation factor related protein 2, ARFRP2, ARF-related protein 2, FLJ20051.

Source:E.coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MSDLRITEAF LYMDYLCFRA
LCCKGPPPAR PEYDLVCIGL TGSGKTSLLS KLCSESPDNV VSTTGFSIKA VPFQNAILNV
KELGGADNIR KYWSRYQGS QGVIFVLSA SSEDDLEAAR NELHSALQHP QLCTLPFLIL
ANHQDKPAAR SVQEIKKYFE LEPLARGKRW ILQPCSLDDM DALKDSFSQL INLLEEKDHE
AVRM

Purity:Greater than 85% as determined by SDS-PAGE.

Formulation:

The ARL15 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 0.1M NaCl and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

ARL15 is a member of the ARF family which are essential in eukaryotic vesicular trafficking pathways and have a vital part in the activation of phospholipase D. ARL15 variants influence levels of Acrp30, an adipocyte-derived protein that is extremely genetic and inversely related to the NeoBiolab of type 2 diabetes mellitus and coronary heart disease.

To place an order, please [Click HERE](#).