

SARS Core (340-390 a.a.)

Description: The E.Coli derived 32kDa recombinant protein contains the Nucleocapsid core protein 340-390 amino acids immunodominant regions.

Catalog #:SAPS-261

Purity: SARS Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

For research use only.

Specificity:

SARS Core protein is Immunoreactive with sera of SARS-infected individuals.

Formulation:

50mM Tris-HCl, 60mM NaCl and 50% glycerol.

Stability:

Five year frozen, 6 month at +40C.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

SARS Core antigen is suitable for ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.

Introduction:

SARS Coronavirus is an enveloped virus containing three outer structural proteins, namely the membrane (M), envelope (E), and spike (S) proteins. Spike (S)-glycoprotein of the virus interacts with a cellular receptor and mediates membrane fusion to allow viral entry into susceptible target cells. Accordingly, S-protein plays an important role in virus infection cycle and is the primary target of neutralizing antibodies.

Storage:

SARS Core protein is shipped at ambient temperature. Upon arrival, store at -200C.

To place an order, please [Click HERE](#).